

PHYSICAL SCIENCE

Product Name: AG-00009


UV Spectrum: 419 nm

Hydrodynamic Size: 113.8 nm

PDI: 0.39



Zeta: -23.65 mV

Average Size: 17.83 ± 9.54 nm




TEM Image of NEM-AG-00009


Size distribution by intensity



Histogram



NEM-AG-00009



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.