

PHYSICAL SCIENCE

Product Name: AG-00026


UV Spectrum: 414 nm

Hydrodynamic Size: 166.7 nm



PDI: 0.44

Zeta: -24.44 mV


Average Size: 36.90 ± 14.84 nm




TEM Image of AMR-AG-00026



Histogram



AMR-AG-00026



It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in **DISCLAIMER:** this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.