

PHYSICAL SCIENCE

Product Name: AG-00059


UV Spectrum: 437 nm

Hydrodynamic Size: 112.8 nm



PDI: 0.35

Zeta: -26.32 mV


Average Size: 23.21 ± 5.76 nm




TEM Image of LAV-AG-00059



LAV-AG-00059



Histogram



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.