

PHYSICAL SCIENCE

Product Name: AG-00063


UV Spectrum: 442 nm

Hydrodynamic Size: 155.9 nm

PDI: 0.23



Zeta: -21.5 mv

Average Size: 76.59 ± 33.14 nm




TEM Image of SJW-AG-00063


Size distribution by intensity



Histogram



SJW-AG-00063



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.