

PHYSICAL SCIENCE

Product Name: AG-00068


UV Spectrum: 410 nm

Hydrodynamic Size: 270.7 nm

PDI: 0.38



Zeta: -27.85 mV

Average Size: 21.64 ± 6.91 nm




TEM Image of BTM-AG-00068


Size distribution by intensity



Histogram



BTM-AG-00068



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.