

PHYSICAL SCIENCE

Product Name: AG-00077


UV Spectrum: 441 nm

Hydrodynamic Size: 158.3 nm




PDI: 0.25

Zeta: -26.29 mV

Average Size: 37.45 ± 10.78 nm



TEM Image of PLS-AG-00077



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.