

PHYSICAL SCIENCE

Product Name: AG-00110


UV Spectrum: 437 nm

Hydrodynamic Size: 153.9 nm




PDI: 0.23

Zeta: -26.44 mV

Average Size: 31.91 ± 11.78 nm



TEM Image of PMT-AG-00110



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.