

PHYSICAL SCIENCE

Product Name: AG-00112


UV Spectrum: 434 nm

Hydrodynamic Size: 117.1 nm




PDI: 0.22

Zeta: -25.32 mv


Average Size: 32.78 ± 16.76 nm



TEM Image of ARK-AG-00112



ARK-AG-00112



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.