

PHYSICAL SCIENCE

Product Name: AG-00116


UV Spectrum: 412 nm

Hydrodynamic Size: 87.76 nm



PDI: 0.48

Zeta: -30.24 mV


Average Size: 24.59 ± 15.02 nm




TEM Image of KRJ-AG-00116



Histogram



KRJ-AG-00116



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.