

PHYSICAL SCIENCE

Product Name: AG-00124


UV Spectrum: 431 nm

Hydrodynamic Size: 386.1 nm

PDI: 0.25



Zeta: -12.61 mV

Average Size: 34.41 ± 26.43 nm




TEM Image of TMR-AG-00124


Size distribution by intensity



Histogram



TMR-AG-00124



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.