

PHYSICAL SCIENCE

Product Name: AG-00127


UV Spectrum: 425 nm

Hydrodynamic Size: 144.6 nm

PDI: 0.33



Zeta: -25.93 mV

Average Size: 34.07 ± 14.06 nm




TEM Image of MTA-AG-00127


Size distribution by intensity



Histogram



MTA-AG-00127



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.