

PHYSICAL SCIENCE

Product Name: AG-00140


UV Spectrum: 429 nm

Hydrodynamic Size: 270.3 nm

PDI: 0.39



Zeta: -9.41 mV

Average Size: 28.88 ± 9.1 nm




TEM Image of ANT-AG-00140


Size distribution by intensity



Histogram



ANT-AG-00140



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.