

PHYSICAL SCIENCE

Product Name: AG-00148


UV Spectrum: 411 nm

Hydrodynamic Size: 56.38 nm

PDI: 0.35


Zeta: -40.25 mV

Average Size: 7 ± 3.24 nm




TEM Image of SKA-AG-00148


Histogram



Size distribution by intensity



SKA-AG-00148



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.