

PHYSICAL SCIENCE

Product Name: AG-00191


UV Spectrum: 430 nm

Hydrodynamic Size: 279.4 nm

PDI: 0.33



Zeta: -27.07 mV

Average Size: 16.13 ± 3.94 nm




TEM Image of BLB-AG-00191


Size distribution by intensity



Histogram



BLB-AG-00191



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.