

PHYSICAL SCIENCE

Product Name: AG-00202


UV Spectrum: 431 nm

Hydrodynamic Size: 79.96 nm



PDI: 0.38

Zeta: -20.6 mV


Average Size: 37.63 ± 16.68 nm




TEM Image of SRM-AG-00202



SRM-AG-00202



Histogram



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.