

PHYSICAL SCIENCE

Product Name: AG-00212


UV Spectrum: 447 nm

Hydrodynamic Size: 267.7 nm

PDI: 0.26



Zeta: -19.05 mV

Average Size: 65.86 ± 35.89 nm




TEM Image of BKB-AG-00212


Size distribution by intensity



Histogram



BKB-AG-00212



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.