

PHYSICAL SCIENCE

Product Name: AG-00225


UV Spectrum: 430 nm

Hydrodynamic Size: 124.3 nm

PDI: 0.29



Zeta: -34 mV

Average Size: 23.19 ± 12.01 nm




TEM Image of ATB-AG-00225


Size distribution by intensity



Histogram



ATB-AG-00225



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.