

PHYSICAL SCIENCE

Product Name: AG-00226


UV Spectrum: 439 nm

Hydrodynamic Size: 655.6 nm

PDI: 0.59



Zeta: -22.7 mV

Average Size: 32.32 ± 18.31 nm




TEM Image of AJW-AG-00226


Size distribution by intensity



Histogram



AJW-AG-00226



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.