

PHYSICAL SCIENCE

Product Name: AG-00269


UV Spectrum: 422 nm

Hydrodynamic Size: 142.6 nm


PDI: 0.25

Zeta: -20.17 mV



Average Size: 28.53 ± 12.33 nm




TEM Image of BAN-AG-00269



Size distribution by intensity



BAN-AG-00269



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.