

PHYSICAL SCIENCE

Product Name: AU-00002

UV Spectrum: 542 nm

Hydrodynamic Size: 223.8 nm



PDI: 0.11

Zeta: -22.29 mV



Average Size: 57.31 ± 16.49 nm



Size distribution by intensity




KDM-AU-00002



TEM Image of KDM-AU-00002

Histogram



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.