

PHYSICAL SCIENCE

Product Name: AU-00031


UV Spectrum: 559 nm

Hydrodynamic Size: 582 nm

PDI: 0.5



Zeta: -14.9 mV

Average Size: 20 +/- 15.9 nm




TEM Image of PPL-AU-00031


Size distribution by intensity



Histogram



PPL-AU-00031



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.