

PHYSICAL SCIENCE

Product Name: AU-00037


UV Spectrum: 544 nm

Hydrodynamic Size: 110.9 nm

PDI: 0.3



Zeta: -11.0 mV

Average Size: 21 +/- 18.9 nm




TEM Image of STV-AU-00037


Size distribution by intensity



Histogram



STV-AU-00037



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.