

PHYSICAL SCIENCE

Product Name: AU-00040


UV Spectrum: 535 nm

Hydrodynamic Size: 77.78 nm

PDI: 0.34


Zeta: -16.56 mV

Average Size: 25.6+/-8.4 nm




TEM Image of JTM-AU-00040


Size distribution by intensity




Zeta potential distribution



JTM-AU-00040



Histogram



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.