

PHYSICAL SCIENCE

Product Name: AU-00041


UV Spectrum: 543 nm

Hydrodynamic Size: 241.4 nm

PDI: 0.23



Zeta: -18.81 mV

Average Size: 13.2 +/- 7.6 nm




TEM Image of SNK-AU-00041


Size distribution by intensity



Histogram



SNK-AU-00041



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.