

PHYSICAL SCIENCE

Product Name: AU-00048


UV Spectrum: 533 nm

Hydrodynamic Size: 68.9 nm



PDI: 0.5

Zeta: -16.4 mV


Average Size: 8.97 +/-6.3 nm




TEM Image of BAL-AU-00048



BAL-AU-00048



Histogram



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.