

PHYSICAL SCIENCE

Product Name: AU-00050


UV Spectrum: 536 nm

Hydrodynamic Size: 118.5 nm

PDI: 0.24



Zeta: -18.4 mV

Average Size: 29.44 ± 20.92 nm




TEM Image of BLV-AU-00050


Size distribution by intensity



Histogram



BLV-AU-00050



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.