

PHYSICAL SCIENCE

Product Name: AU-00052


UV Spectrum: 543 nm

Hydrodynamic Size: 190.6 nm

PDI: 0.20



Zeta: -16.91 mV

Average Size: 239.02 ± 18.53 nm




TEM Image of WTG-AU-00052


Size distribution by intensity



Histogram



WTG-AU-00052



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.