

PHYSICAL SCIENCE

Product Name: AU-00055


UV Spectrum: 541 nm

Hydrodynamic Size: 352.6 nm

PDI: 0.36



Zeta: -21.78 mV

Average Size: 16.86 ± 6.67 nm




TEM Image of FEN-AU-00055


Size distribution by intensity



Histogram



FEN-AU-00055



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.