

PHYSICAL SCIENCE

Product Name: AU-00064


UV Spectrum: 549 nm

Hydrodynamic Size: 243.9 nm

PDI: 0.21



Zeta: -20.76 mv

Average Size: 24.70 ± 10.71 nm




TEM Image of MIT-AU-00064


Size distribution by intensity



Histogram



MIT-AU-00064



It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in **DISCLAIMER:** this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.