

PHYSICAL SCIENCE

Product Name: AU-00065


UV Spectrum: 542 nm

Hydrodynamic Size: 329.9 nm

PDI: 0.57



Zeta: -8.63 mv

Average Size: 34.39 ± 18.28 nm




TEM Image of MSW-AU-00065


Size distribution by intensity



Histogram



MSW-AU-00065



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.