

PHYSICAL SCIENCE

Product Name: AU-00069


UV Spectrum: 534 nm

Hydrodynamic Size: 84.51 nm




PDI: 0.48

Zeta: -12.52 mV

Average Size: 71.17 ± 35.08 nm



TEM Image of BLB-AU-00069



It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in **DISCLAIMER:** this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.