

PHYSICAL SCIENCE

Product Name: AU-00078


UV Spectrum: 554 nm

Hydrodynamic Size: 130.4 nm

PDI: 0.27



Zeta: -20.03 mv

Average Size: 38.04 ± 18.72 nm




TEM Image of PSP-AU-00078


Size distribution by intensity



Histogram



PSP-AU-00078



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.