

PHYSICAL SCIENCE

Product Name: AU-00092


UV Spectrum: 531 nm

Hydrodynamic Size: 179.6 nm

PDI: 0.23



Zeta: -7.7 mV

Average Size: 16.90 ± 5.88 nm




TEM Image of SFD-AU-00092


Size distribution by intensity



Histogram



SFD-AU-00092



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.