

PHYSICAL SCIENCE

Product Name: AU-00101


UV Spectrum: 538 nm

Hydrodynamic Size: 215.3 nm

PDI: 0.47



Zeta: -6.62 mV

Average Size: 36.26 ± 10.13 nm




TEM Image of CRA-AU-00101


Size distribution by intensity



Histogram



CRA-AU-00101



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.