

PHYSICAL SCIENCE

Product Name: AU-00107


UV Spectrum: 532 nm

Hydrodynamic Size: 120.6 nm

PDI: 0.48



Zeta: -22.43 mV

Average Size: 29.85 ± 3.81 nm




TEM Image of CME-AU-00107


Size distribution by intensity



Histogram



CME-AU-00107



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.