

PHYSICAL SCIENCE

Product Name: AU-00129


UV Spectrum: 539 nm

Hydrodynamic Size: 90.08 nm

PDI: 0.27



Zeta: -13.69 mV

Average Size: 20.69 ± 8.32 nm




TEM Image of BLP-AU-00129


Size distribution by intensity



Histogram



BLP-AU-00129



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.