

PHYSICAL SCIENCE

Product Name: AU-00130


UV Spectrum: 536 nm

Hydrodynamic Size: 87.34 nm




PDI: 0.28

Zeta: -21.28 mv


Average Size: 23.97 ± 11.73 nm



TEM Image of KGH-AU-00130



KGH-AU-00130



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.