

PHYSICAL SCIENCE

Product Name: AU-00131


UV Spectrum: 541 nm

Hydrodynamic Size: 171.1 nm

PDI: 0.23



Zeta: -14.08 mV

Average Size: 16.91 ± 7.31 nm




TEM Image of STA-AU-00131


Size distribution by intensity



Histogram



STA-AU-00131



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.