

PHYSICAL SCIENCE

Product Name: AU-00134


UV Spectrum: 534 nm

Hydrodynamic Size: 198.1 nm

PDI: 0.18



Zeta: -29.46 mV

Average Size: 30.41 ± 10.27 nm




TEM Image of PGF-AU-00134


Size distribution by intensity



Histogram



PGF-AU-00134



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.