

PHYSICAL SCIENCE

Product Name: AU-00141


UV Spectrum: 534 nm

Hydrodynamic Size: 39.06 nm



PDI: 0.51

Zeta: -23.34 mV


Average Size: 30.78 ± 11.84 nm



TEM Image of MNL-AU-00141



MNL-AU-00141



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.