

PHYSICAL SCIENCE

Product Name: AU-00178


UV Spectrum: 558 nm

Hydrodynamic Size: 243 nm

PDI: 0.22



Zeta: -17.39 mV

Average Size: 22.22 ± 8.29 nm




TEM Image of LKP-AU-00178


Size distribution by intensity



Histogram



LKP-AU-00178



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.