

PHYSICAL SCIENCE

Product Name: AU-00183


UV Spectrum: 559 nm

Hydrodynamic Size: 99.71



PDI: 0.44

Zeta: -29.34 mV


Average Size: 48.69 ± 26.08 nm




TEM Image of SLK-AU-00183



SLK-AU-00183



Histogram



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.