

PHYSICAL SCIENCE

Product Name: AU-00192


UV Spectrum: 531 nm

Hydrodynamic Size: 555 nm

PDI: 0.75



Zeta: -26.91 mv

Average Size: 35.99 ± 17.38 nm




TEM Image of HEN-AU-00192


Size distribution by intensity



Histogram



HEN-AU-00192



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.