

PHYSICAL SCIENCE

Product Name: AU-00209


UV Spectrum: 537 nm

Hydrodynamic Size: 114 nm



PDI: 0.44

Zeta: -20.05 mV


Average Size: 32.95 ± 14.96 nm




TEM Image of CNM-AU-00209



Histogram



CNM-AU-00209



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.