

PHYSICAL SCIENCE

Product Name: AU-00211


UV Spectrum: 543 nm

Hydrodynamic Size: 88.28 nm


PDI: 0.3177

Zeta: -28.31 mV


Average Size: 18.37 ± 6.09 nm



TEM Image of MGW-AU-00211



Histogram



MGW-AU-00211



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.