

PHYSICAL SCIENCE

Product Name: AU-00213


UV Spectrum: 531 nm

Hydrodynamic Size: 79.91 nm

PDI: 0.47



Zeta: -13.89 mV

Average Size: 34.24 ± 15.63 nm




TEM Image of BKB-AU-00213


Size distribution by intensity



Histogram



BKB-AU-00213



DISCLAIMER: It is common for metallic nanoparticles to agglomerate with time. Therefore, sometimes, the sizes as reported in this catalogue may not be reproducible upon purchase. We strongly recommend sonication before use to minimize agglomeration.